Appareils de laboratoire

Bain-marie - étuves

Laboratory equipment

Centrifuges - ovens


**Ref :
591 004**

Français – p 1

English – p 11

Version : 6010

Bain marie à cuve transparente
Water bath with transparent tank

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
NOTES

The device is designed for experiments that require a preparation to be maintained for different intervals at a given temperature, between the ambient temperature and 80°C.

It is designed to hold up to 12 (180 x 18 mm) test tubes (minimum height of tubes: 130 mm).

1 Description

The device consists of three removable parts:

1.1 The heating resistance and the control box

The heating resistance is in the form of a M-shaped folded tubing in its horizontal heating portion. This ensures homogeneous heating. Parallely, a tube with a closed end, containing the temperature probe (thermistor), is installed on the vertical portion of this tubing, the non-heating part that corresponds to the power supply of the resistance.

The control box, made of resistant plastic (Dimensions: 140 x 95 x 30 mm) includes:

- On the front face: The "On-Off" switch and the temperature adjustment button
- On top, A POWER-ON RED INDICATOR LAMP and A GREEN INDICATOR LAMP that comes on only when the resistance is heating, and goes off when the displayed temperature is reached
- On the lower portion of the box: the mains cord 2P + T.

The heating resistance and the control box are connected together by a plastic block that fits on one of the sides of the tank. The geometry of the assembly makes it extremely stable.

The device is equipped with a float-based minimum water level detector that prevents heating when there is no water. In this case, the red indicator lamp comes on but, whatever be the settings, the resistance does not heat and the green indicator lamp does not come on.

1.2 The transparent tank

It contains the liquid to be heated.

The minimum water level is screen-printed on one of its main faces.

Its inner dimensions are: 105 x 157 x 140 mm.

A block with the control box and the resistance fits onto one of its smaller sides.

1.3 The test tube support

It can hold 12 test tubes of \varnothing 18 mm maximum and a control thermometer (not supplied with the device).

2 Technical specifications

- Supply voltage: 215 to 230 Volts with ground terminal
- Power: 600 W at 220 Volts
- Capacity of the tank: About 2 litres
- Overall dimensions: 215 x 110 x 160 mm
- Mass : About 1.2 kg
- Temperature regulation: from 20 to 65°C: $\pm 1^\circ\text{C}$ then 1.5°C up to 80°C
- Reading accuracy : from 20 to 65°C: $\pm 2.5^\circ\text{C}$ then 3°C up to 80°C

2.1 Electronic regulation

The water bath has an electronic temperature control. The electronic regulator eliminates the production of radio interferences and mains distortions.

3 Operation – commissioning

3.1 Implementation of the device

1°) For safety reasons, make sure that the socket to which the device will be connected has a proper ground terminal.

2°) Position the adjustment button on the desired temperature setting.

3°) Connect the water bath to the mains (2P + T). It has been designed to operate at 220 Volts. At this voltage, the heating capacity is correct. Powered at 110 Volts, the heating capacity reduces by up to 4 times, thus it takes much longer to attain the desired temperature. However, the electronic regulation and the heating resistance control operate correctly up to 75°C.

4°) Immerse the heating resistance into the tank by making sure that there is sufficient water (demineralised water to avoid calcification on the tank and the resistance).

5°) Put the mains switch "On".

The red indicator lamp comes on.

The green indicator lamp comes on if the resistance is heating. There are three cases:

- The green indicator lamp is on: the heating is constant, the temperature is less than that displayed
- It comes on and goes off periodically: the temperature is quite close to the desired value
- The indicator lamp goes off: the bath is at the desired temperature (or slightly above).

6°) If the water level is not sufficient, the "lack of water" safety system prevents heating and the green indicator lamp does not come on.

PRECAUTIONS FOR USE

All metal parts of the device are connected to the "Ground" for safety reasons.

***The control box is not watertight.
As a result:***

NEVER IMMERSE THE CORD OR THE CONTROL BOX IN ANY LIQUID

The device must mandatorily be connected to a network that has a ground fault circuit interrupter. If this is not available, install an isolating transformer.

3.2 Controlling the operation and fine adjustment

The figure below gives the curves for increase of temperature for 3 positions of the thermostat:

These show that we enter the regulation zone in a time that varies between 5 and 15 minutes depending on the temperature displayed, when $V = 220$ Volts and when there are 1.6 litres of water in the tank.


It must be observed that for relatively low temperatures, the indicator lamp going off does not indicate that the regulation is being done correctly: release of the heat stored in the resistance body makes the temperature of the bath rise further: exchanges with the outside being that much slower as the temperature is close to the ambient temperature.

This is why, when the bath needs to be maintained at a very precise temperature, it is advisable to first adjust the setting button to a temperature that is slightly below the desired temperature (3 to 4°C). "Additional" corrections may be made when the regulation is well established.

As in any regulated device, there may be a slight difference between the temperature displayed and the temperature read on the thermometer inside. This difference may be due to two main reasons: the non-linearity of the display potentiometer and above all the temperature gradient of the liquid in the bath.

If we need to perform repetitive operations or if we require a more precise adjustment, we will note once and for all the corresponding position of the button or draw a nomogram indicating the corrections to be made.

When we know our device well, this correction becomes rapidly automatic for the operator.

4 Care – maintenance – repairs

4.1 Maintenance of the tank

- It is recommended that demineralised water be used in the tank to avoid any calcification on the tank itself, on the thermometer and on the heating resistance.
- To clean the tank, use only water, without any other product (alcohol, acetone, trichloroethylene corrode matter). A weak acid may however be used.
- Do not scrape the tank, not even with a rounded blade or an abrasive buff (risk of scratches).

4.2 Maintenance of the resistance

Use ordinary detergent to clean it.

4.3 In case of failure to operate

If after having connected the device to the mains and switched it on, we do not observe any heating:

1°) Make sure that the red lamp is on. Otherwise, check that there is voltage at the mains socket.

2°) Check that the thermostat does not indicate a temperature that is less than that of the water bath.

If for any reason the case needs to be opened, disconnect the mains cord.

5 Experiments

One must first remember that the device must never be left powered on when the resistance is not immersed in water.

Make sure that the water reaches at least the minimum level indicated on the tank.

If the heating resistance needs to be taken out of the water while it has still not cooled down, do not touch the heating element with the fingers, or place it on a surface that is sensitive to heat.

5.1 Examples of use of the water bath

Uses of a water bath in Biology are too many for all of them to be mentioned.

We will limit ourselves here to suggesting some of the standard experiments that may interest secondary school classes. The most common ones have been underlined.

6 Study of carbohydrates

6.1 Reduction of silver salts by glucose

Temperature: 50 to 60° .

The experiment will be performed on a silver ammonium nitrate solution obtained by adding a quantity of pure diluted ammonia water, just enough to remove the precipitate that forms at the start (AgOH precipitate), to a 2% silver nitrate solution.

Then, a few ml of pure glucose solution at 0.2% will be added and the mixture will be placed in the water bath for a few minutes. To obtain homogenous silvering that adheres properly, it is necessary to use a test tube that is perfectly clean (rinsed, if necessary, with dilute nitric acid, then with alcohol and finally with distilled water) and pure products.

Under the same conditions, saccharose does not give any result: saccharose is a non-reducing sugar.

6.2 Hydrolysis of starch by ptyalin (saliva)

The water bath 591004 is used to process 12 test tubes simultaneously and to follow the progress of the reaction over time.

- a) Prepare starch paste to be added to 100ml of water brought to the boil, a suspension of 1 g of potato starch in 10 ml of water. Let it boil for 1 to 2 minutes.
- b) Place about 5 ml of starch paste in each test tube and add about 1 ml of saliva.
- c) Test the first tube with Lugol's solution: blue coloration shows the presence of unconverted starch.
- d) The 11 other test tubes are placed in the water bath at 40° and every 10 minutes one tube is removed and tested with iodised water. We obtain a range of colours.
- e) In the last tube, we add Fehling's solution and bring it to boil on a Bunsen burner. The red colouring obtained shows the formation of a reducing disaccharide. We verify that only saliva does not reduce Fehling's solution.

6.3 Effect of temperature on enzyme action

2 water baths, which include one at 80°C, are required so that the experiment does not last too long.

Prepare 16 test tubes with 1 ml of distilled water and 2 drops of Lugol in each of them.

1 tube will be marked "0°C", another "+ 80°C".

Six tubes will be marked "ambient" and six tubes "50°C".

2 ml of starch paste will be placed in each of the tubes.

The 0°C tube will be placed in melting ice, while the 20° tubes will be placed in a tube holder. The two other sets will be placed in the two water baths set at 50°C and 80°C respectively.

Add 1 ml of saliva to the "0°C" tube and test immediately and after half an hour.

Do the same with the other tubes (only the "100°C" tube to be tested every 5 minutes).

Starch paste and saliva that has been brought to boil beforehand will be added to the 80°C tube and this will be tested immediately and after half an hour. The colouring ranges obtained will be compared and grouped together by similar colouring.

We will thus observe that we obtain identical colourings for tubes processed at different temperatures but also for different time intervals. The results must be represented in the form of a graph.

We could take the 0° tube once again, place it at 40° for example, and test it after 15 minutes. Lugol will show that hydrolysis has occurred: at 0° the effect of ptyalin is suspended but in a reversible manner.

On the other hand, the 80°C tube, even when brought to ambient temperature, will not show any enzymatic action: the enzyme has been destroyed by heat.

7 Study of lipids

Saponification of milk cream by pancreatin

Place 10 ml of untreated milk, 1 ml of litmus solution and a drop of soda lye in two test tubes.

Add 3 ml of pancreatin in one of the tubes, and 3 ml of pancreatin that has been boiled beforehand in the other tube.

Place both tubes in the water bath at 40°C for 30 min: the first tube turns pink due to the appearance of fatty acids, while there is no change in the other tube (inhibition of pancreatin due to heat).

8 Study of proteins

8.1 Action of pepsin on fibrin or casein:

Fibrin may be recovered by washing the clot after a lesson on blood. The experiment can also be done with gelatine.

Place the same quantity of proteins in four test tubes. Add:

- In tube 1: 5 ml of pepsin at 0.5 % and 5 ml of HCl at 0.4 %
- In tube 2: 5 ml of pepsin and 5 ml of water
- In tube 3: 5 ml of pepsin boiled beforehand and 5 ml of HCl at 0.4 %
- In tube 4: 5 ml of water and 5 ml of HCl at 0.4 %.

Place all four tubes in the water bath at 40°C for 40 minutes. At the end of this period, only the contents of tube 1 have undergone a significant transformation and the biuret test shows that pepsin has caused a part of the protein to be transformed into polypeptides.

8.2 Action of pH on the hydrolysis of pancreatin

Put equal quantities of fibrin (for example) into four test tubes.

Add:

- Tube no. 1: 3 ml of pancreatin solution at 1 % and 3 ml of water.
- Tube no. 2: 3 ml of pancreatin solution at 1 % and a drop of phenolphthalein and a sodium carbonate solution at 0.5 % until an incipient pink colour is obtained.
- Tube no. 3: Same operation by increasing the amount of carbonate until a clear pink colour is obtained.
- Tube no. 4: 3 ml of pancreatin solution and 3 ml of boric acid solution at 2 %.

Place all four tubes in the water bath at 40°C for 40 minutes. The hydrolysis is complete only in tube 1 (optimum pH) releasing amino acids that neutralise the carbonate (disappearance of the pink colour).

9 Study of haemoglobin

Release of oxygen fixed on oxyhaemoglobin

(Caution: the water bath support is not suitable for this experiment).

Potassium ferricyanide transforms oxyhaemoglobin into methaemoglobin incapable of fixing oxygen which, when liberated, comes out in the form of bubbles.

Pour 30 ml of defibrinated blood into a 180 x 18 test tube and add 1 ml of distilled water to it.

Place 10 ml of a saturated solution of potassium ferricyanide into an ordinary test tube (160 x 16).

Place the two test tubes for 10 minutes in the water bath at 50°C, then pour the ferric solution into the large tube. We can clearly see oxygen bubbles along the walls when the tube is tilted.

10 Study of catabolism

10.1 Effects of temperature on alcoholic fermentation and determination of the products obtained

Dissolve 15 g of glucose in 50 ml of water and dilute in it about 0.5 g of baker's yeast.

Place the mixture in a 100 ml conical flask (Erlenmeyer) stopped with a cork having a hole fitted with a pressure relief tube with 1 bend.

Fill up to half a 100 or 200 ml capacity flask (straight neck or conical) with limewater. Plug the container with a cork having two holes, one of them fitted with a pressure relief tube with 1 bend.

Place the vial containing the sugar solution in the water bath, the other flask remaining on the table.

Connect the two pressure relief tubes with a rubber or flexible plastic tubing.

The reaction begins rapidly and when the gas output seems constant, count the number of bubbles emitted in a given time (one minute, for example).

Power on the water bath and set to 30°. After stabilisation, count once again the number of bubbles emitted and continue in this manner by increasing the temperature in plateaus of 10° up to 80°).

Draw a graph by plotting temperature on the x-axis and the number of bubbles on the y-axis. There is a maximum at about 40° and at 60° there is no emission.

The turbidity of the limewater demonstrates the emission of CO₂ and the presence of alcohol in the liquid is detected by green colouring if a measure of potassium dichromate and a few drops of sulphuric acid are added to it.

10.2 Demonstrating a dehydrogenating enzyme in living tissues.

This is an oxidation-reduction experiment, which shows that there exists an enzyme capable of causing an oxidation by dehydrogenation in living tissues. The indicator will be a body capable of fixing the hydrogen released: this is the case with methyl blue that gets discoloured.

Fill 5 (180 x 18 mm) test tubes half with Ringer's solution:

- NaCl : 6.5 g
- KCl : 0.14 g
- $\text{PO}_4 \text{NaH}_2$: 0.20 g
- CaCl_2 : 0.12 g
- CO_3NaH : 0.20 g
- Distilled water : 1,000 g

Introduce the body of two mussels (without the shell) in tubes 1 and 2.

Place a piece of turnip cut with a cork borer in tubes 3 and 4.

Place tubes 2 and 4 in boiling water for 5 minutes.

Add to each tube 5 drops of methyl blue at 0.1 % then paraffin oil that will form a screen in relation to the air. Tube 5 without living tissue will be used as control sample.


Place the five tubes in the water bath at 35 ° for half an hour and compare the colours: only tubes 1 and 3 show a discoloration. The deoxydase present in the tissues of the mussel and the turnip is destroyed by heat.

11 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

Laboratory equipment
Water bath with transparent tank
Ref :
591 004

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service * :
+ 33 (0)2 32 29 40 50

** Call cost not included*