Méthode VTT[®]

VTT[®] method

**Ref :
472 050**

Français – p 1

English – p 5

Version : 6005

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VTT[®] Charge Sensor

1 GENERALITES

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1 GENERAL INFORMATION

This sensor is used to measure the electrostatic charge with the VTT® system.

The sensor is supplied with a shielded cord equipped with 2 alligator clips.

2 USE

2.1 Interfaces

This sensor works with the VTT® console or ESAO Primo.

2.2 Connecting the sensor to the VTT® console

The sensor is connected using the connector located at the back of the box of the VTT console or the front face of the ESAO Primo:

- Place the sensor in line with the selected measurement channel.
- Engage the sensor on the console guide.
- Follow the guidance system of the VTT® console.
- Engage it fully.

2.3 Commissioning with the VTT console

- Connect the sensor.
- Switch on the VTT® console.

If the sensor is not detected, update your VTT console connected to the PC using the floppy disc via the Espace VTT software.

- Connect the console to the PC,
- Open the Espace VTT by clicking on the icon  present in "My computer",
- Click on Install ,
- Install the Biosflash software.

2.4 Commissioning with the ESAO Primo

The sensor is automatically detected by the sensor.

- Connect the console to your PC,
- Connect the sensor to the console,
- Start your acquisition software.

2.5 Scale

The user has a measurement scale of -100 nC / + 100 nC. (nano-coulomb)

3 EXAMPLES OF EXPERIMENTS

3.1 Electrification by contact:

1) Discharge the sensor by pressing the Reset button.

Rub an ebonite rod with cabrettas. The rod gets negatively charged. Bring the end of the rod close to the red clip until contact. A part of the negative charges is transferred. The sensor gets negatively charged. The (-) sign and the value of the charge in nC remain displayed even when the rod is taken away.

2) Discharge the sensor by pressing the Reset button.

Rub an ebonite rod with cabrettas. The rod gets negatively charged. Place the end of the rod, then the entire length of the rod in contact with the red clip.

Compare the 2 values.

(To predict the sign of the charge, refer to the triboelectricity list on the next page)

3.2 Electrostatic induction:

1) Discharge the sensor by pressing the Reset button.

Rub an ebonite rod with cabrettas. Bring the rod close to the red clip without touching it. The (-) sign and the value of the charge in nC are displayed only when the rod is close to the clip.

2) Discharge the sensor by pressing the Reset button.

Rub an ebonite rod with cabrettas. Bring the rod close to the red clip and press the reset button. The free electrons of the clip are repelled by the negative charges of the rod and discharged through the ground. By removing the rod and releasing the Reset button, the sensor remains charged positively because it has lost electrons. The (+) sign and the value of the charge in nC are displayed.

Triboelectric series

	Acetate
	Rabbit fur
	Glass
	Dry human hair
	Nylon
MORE POSITIVE (+)	Wool
	Cat fur
	Cabrettas
	Lead
	Silk
	Aluminium
	Paper
	Cotton
	Wood
	Amber
	Steel
	Ebonite
	Nickel - Copper
MORE NEGATIVE (-)	Brass – Stainless steel
	Rubber
	Most thermoplastics (PVC, etc.)
	PTFE (polytetrafluoroethylene, Teflon)

When we rub a material **A** with material **B** placed below it in the list, **A** gets charged with positive electricity and **B** gets charged with negative electricity.

For example, glass rubbed with wool, cabrettas, cotton or nylon will always be positive. If it is rubbed on rabbit fur, then the glass will be negative.

At the top of the list we find materials that have a greater tendency to shed electrons and at the bottom, materials that are most "avid" for electrons. The intensity of a triboelectric charge depends on the following factors:

- Separation between the two materials in the triboelectric table.
- Experimental conditions (temperature, relative humidity).
- Conditions of the material surfaces (cleanliness, texture).
- Contact pressure during rubbing.
- Speed of separation of materials after rubbing.

4 TECHNICAL SPECIFICATIONS

- Measuring range: from -100 nC to + 100 nC
- Resolution: 0.2 nC
- Maximum operating voltage: 30 V

5 AFTER-SALES SERVICE

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts, please contact:

JEULIN – TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries