TP de chimie

Chemistry tutorial

**Ref :
282 072**

Français – p 1

English – p 7

Version : 6010

Ensemble étude de piles

Cells study kit

1 Description

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
1 Description

This complete kit allows studying up to 10 cells by combining 5 oxidation-reduction pairs (Ag^+/Ag , Cu^{2+}/Cu , Pb^{2+}/Pb , Fe^{2+}/Fe , Zn^{2+}/Zn). The measurements of intensity and voltage at the terminals of the cell are obtained using a multimeter.

The support keeps the different components in place for a greater comfort in usage: the students concentrate only on the measurements to be obtained.

- | | |
|---|---|
| 1 – Support made of plastic material. | 5 – Filter paper sheet. |
| 2 – Slots for the 30 mL flasks and 25 mL beakers | 6 – Water-proof label. |
| 3 – Banana safety plugs $\varnothing$ 4 mm. | 7 – Glass flasks of 30 mL |
| 4 – Electrode holder with rear outlet mounted on banana safety plug with rear outlet $\varnothing$ 4 mm | 8 – Set of electrodes (zinc, iron, lead, copper, silver) identified by a colour code. |

2 Description

This apparatus allows studying the cells by combining the 5 oxidation-reduction pairs available. The set up is made on a plate made of plastic material holding the accessories in place for a more comfortable usage.

The plate is designed for accepting, in its two slots, flasks of 30 mL (delivered with the kit) or beakers of 25 mL (to be ordered separately). The saline bridge between the half-cells is created with a strip of filter paper that provides the electrical conduction, either by soaking it in an ionic solution (potassium nitrate), or by allowing the two solutions of flask to be joined by capillarity. Working with small volumes allows economising the metallic salt solutions that can be used several times.

2.1 Composition:

- A support in plastic material resistant to corrosion
- Two electrode-holders equipped with banana safety plugs with rear outlet for connecting a multimeter
- A set of six metallic electrodes (iron, copper, zinc, lead, silver, tin⁽¹⁾)
- A set of 30 mL flasks for containing the solutions of metallic salts
- A sheet of filter paper 400 × 500 mm to be cut into strips
- A sheet of water-proof labels for identifying the flasks and the stoppers.
- Detailed pedagogical notes.

The metallic electrodes are identified by a colour code, as described in the following table (link between colour and nature of the electrode):

Metal	Colour
Zinc	Black
Iron	White
Tin	Green
Lead	Red
Copper	-----
Silver	Violet

Table 2

The kit is delivered with all the accessories necessary for the experiment except the decimolar solutions of metallic salts.

The decimolar solutions (not supplied) should be prepared from metallic salts corresponding to the nature of electrodes used. For selecting the metallic salts suitably, refer to the instructions in the following chapters.

(1) Due to the toxicity of lead salts, an electrode in tin is supplied for allowing the user to study the pair Sn^{2+}/Sn instead of the pair Pb^{2+}/Pb (refer Chapter 4).

3 Accessories and consumables

Multimeter CL 3005	Ref. 291 156
Beaker 25 mL low shape (set of 10)	Ref. 713 502
Ammonium iron (II) sulphate hexahydrate pure 250 g	Ref. 107 259
Silver nitrate pure 20 g	Ref. 104 030
Copper (II) sulphate pentahydrate pure 250 g	Ref. 105 028
Lead (II) nitrate pure 250 g	Ref. 104 018
Potassium nitrate pure 500 g	Ref. 104 010
Zinc sulphate pure monohydrate 250 g	Ref. 105 142
Demineralised water 1L	Ref. 107 340

4 Usage instructions and precautions

4.1 The silver salts

The solutions of silver salts being photosensitive, they should be kept in brown glass flasks.

It is advisable not to use solutions of metallic chlorides for preventing the solution of silver nitrate via the saline bridge (formation of a precipitate of silver chloride). Use rather solutions of metallic sulphates and nitrates.

4.2 The iron (II) salts

The solutions of iron (II) salts have a limited life period due to the oxidation of iron (II) ions into iron (III) in contact with oxygen in the air. Therefore, it is preferable to use a solution of Mohr salts (ammonical iron (II) sulphate, $(\text{NH}_4)_2\text{Fe}(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}$) whose stability is better than that of iron (II) salts such as the ferrous sulphate (II) or the ferrous nitrate (II). If, however, the colour of the solution turns from light green to brown – characteristic colour of iron (III) ions – it is necessary to replace with a freshly prepared solution.

4.3 The lead salts ⁽¹⁾

Their usage presents a danger for the health of the user and the environment.

T toxic

N dangerous for the environment

Gloves

Eye glasses

Mask

It is necessary to use all the personal and collective protection equipment (gloves, eye glasses...) for preventing any contact by inhalation ⁽²⁾, by contact with the skin and by ingestion. It is also necessary to prevent the contamination of the working environment and the residues must not be dumped into the environment (sewer, ground..).

Before using these products, consult their Safety Data Sheets on our Internet site:

www.jeulin.fr

(1) The oxidation-reduction potential of the couple Pb^{2+}/Pb being close to that of the couple Sn^{2+}/Sn (reference: ESH $E^\circ \text{Sn}^{2+}/\text{Sn} = -0,14 \text{ V}$, $E^\circ \text{Pb}^{2+}/\text{Pb} = -0,13 \text{ V}$), it is possible to replace the latter with tin (II) salts (SnSO_4 , for example). These are far less toxic. They are however less soluble in the water with neutral pH and the shelf life of tin salts in the solution is limited.

(2) The use of lead salts in aqueous solution eliminates the risks relating to their inhalation.

5 Procedures and measurements

5.1 The traditional experiment


Its photograph is shown on the right.

A cell is created by the association of two half cells using the oxidation-reduction pairs Zn^{2+}/Zn and Cu^{2+}/Cu .


A saline bridge connects the two flasks containing the ionic solutions and the electrodes.

A digital voltmeter measures the e.m.f. of the cell.


5.2 The theoretical operation of the cell


Traditionally, a cell is represented as follows:

For example, for the zinc/copper cell shown above:


On the copper electrode (terminal +), the following reduction is produced:


In the solution, the Cu^{2+} ions that are reduced are replaced by the K^{+} ions of the ionic bridge.

At the zinc electrode (terminal -), the following oxidation is produced:

In the solution, the appearance of Zn^{2+} ions is compensated by the migration of NO_3^{-} ions from the ion bridge.

The energy exchanges in the cell and in the exterior circuit are as follows:

The positive pole of the cell is associated with the electrodes of the oxidation-reduction pair containing the stronger oxidant: Cu^{2+} .

The negative pole of the cell is associated with the electrode of the oxidation-reduction pair containing the stronger reducer: Zn.

The e.m.f. of the cell studied is $E = 1.08 \text{ V}$.

5.3 The study of other cells


The following table gives the measured values of all the pairs of cells that can be created with this kit:

Voltmeter	The values below are expressed in volt				
Terminal – (Com) Terminal +	Ag	Cu	Pb	Fe	Zn
Ag	0	0.45	1	1.1	1.52
Cu	-0.44	0	0.50	0.73	1.08
Pb	-0.97	-0.51	0	0.26	0.59
Fe	-1.15	-0.72	-0.29	0	0.35
Zn	-1.53	-1.08	-0.62	-0.36	0


(The sign of the voltage measured corresponds to the connection indicated of the voltmeter).

6 Classification of M^{n+}/M pairs

By choosing arbitrarily the value 0 for the pair Cu^{2+}/Cu , the experimental classification of metallic oxidation-reduction pairs is established as follows:

The cell whose e.m.f. is the highest (1.53 V) is given by the following representation:

Note: These experimental values may vary from theoretical values.

7 After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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Service gratuit * :
+ 33 (0)2 32 29 40 50

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Rubrique FAQ

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

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