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Measurement and instrumentation

**Ref :
293 057**

Français – p 1

English – p 17

Version : 9001

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NOTES


- ① On – off switch
- ② Power-on indicator lamp
- ③ Voltage selector switch: $3 \times 6 \text{ V} - 3 \times 9 \text{ V}$ RMS between the live and neutral
- ④ Frequency selector switch: 0.1 Hz - 1 Hz - 50 Hz -150 Hz and “0.1 Hz” with manual displacement
- ⑤ Output terminals

1 Measuring the phase to ground and phase to phase voltages

In a three-phase installation there are three lives and one neutral wire.

The voltage between the phase and the neutral is called phase to ground voltage; its effective value is noted as: V .

The voltage between the two phases is called the line voltage; its effective value is noted as: U .

1.1 Objective

To check the relation between phase to ground voltage V and line voltage U :

$$U = V \times \sqrt{3}.$$

1.2 Equipment

- An extra-low voltage three-phase generator;
- Three digital voltmeters
- Six connecting wires.

1.3 Measurements

- Measure the phase to ground and line voltages; complete the table of measurements.

Example of results:

Rated capacity	Phase to ground voltages			Line voltages		
$3 \times 6 \text{ V}$	V_{1N}	V_{2N}	V_{3N}	U_{12}	U_{23}	U_{31}
$3 \times 9 \text{ V}$	$V_{1N} = 8.9 \text{ V}$	$V_{2N} = 8.9 \text{ V}$	$V_{3N} = 8.9 \text{ V}$	$U_{12} = 15.5 \text{ V}$	$U_{23} = 15.3 \text{ V}$	$U_{31} = 15.4 \text{ V}$

- Compare for each of the two rated capacities the values of the phase to ground voltages V then the line voltages U .
- Check the relation $U = V \times \sqrt{3}$

2 Viewing the voltages using ESAO4 and Générís

2.1 Objective


- To view in three-phase sinusoidal mode: The voltages between the phase and the neutral
The voltages between the phases.
- To verify that a three-phase voltage system is a system with three sinusoidal voltages having the same amplitude, same frequency and offset one from the other by a third of a period
- To verify that the line voltages are leading by $\frac{\pi}{6}$ from the phase to ground voltages
- To construct Fresnel's diagram for the phase to ground and line voltages (see 4).

2.2 Equipment

- An extra-low voltage three-phase generator (9V - 1 Hz)
- A computer with Générïs installed
- ESAO4 console
- Three voltmeter adaptors
- Six connecting wires.

Phase to ground voltages


2.3 Set-up

2.4 Experiment

- Perform the assembly.
- Start the computer; click on the Générïs icon (the indicator lamp on the console comes on).
- Select the Test and calibrate icon
 - o Select: Channel: **Voltmeter /**
 Rated capacity: Ⓢ 25 V
 Ⓢ Instantaneous
 - o Adjust the zero
 - o Click on OK
 - o Select: Channel: **Voltmeter//**
 Rated capacity: Ⓢ 25 V
 Ⓢ Instantaneous
 - o Adjust the zero
 - o Click on OK
 - o Select: Channel: **Voltmeter///**
 Rated capacity: Ⓢ 25 V
 Ⓢ Instantaneous

- o Adjust the zero
- o Click on OK
- Select the Parameters icon
 - o Validate: Acquisition mode: f(t) ▼
 Acquisition time: 2 s ▼
 Number of points: 501 ▼
☒ Synchronisation
 - o Select: Channels on the ordinate: Voltmeter/
 Voltmeter//
 Voltmeter///
 - o Click on Synchro...
 - o Select: Channel: Voltmeter/ ▼
 Level: 0
☒ Synchronisation ☐ Increasing
 - o Click on OK
 - o Click on OK
 - Select the Launch icon
 - o Select: Magnitude: v1
 Names of the curves to be acquired: v1(t)
 Magnitude: v2
 Names of the curves to be acquired: v2(t)
 Magnitude: v3
 Names of the curves to be acquired: v3(t)
 - o Click on Acquisition

Carry out the same search for:

$$v_2(t) = a \cdot \sin(w \cdot t + b)$$


and

$$v_3(t) = a \cdot \sin(w \cdot t + b)$$

We obtain for example:

Parameters a ☒ 12.8
 w ☒ 6.29
 b ☒ 2.02
 Standard deviation S =
 201mV

Parameters a ☒ 12.9
 w ☒ 6.27
 b ☒ 4.14
 Standard deviation S =
 209mV

Check that $v_1(t) + v_2(t) + v_3(t) = 0$ ($\vec{V}_1 + \vec{V}_2 + \vec{V}_3 = \vec{0}$).

- Select Processing then Formula
 - o Select: Curve to be processed: v1(t) ▼

Ordinate:	Name: v
	S.I. Units: V
Formula	v = : v1+v2+v3
Resulting curve:	v(t)
 - o Click on Plot then Save and then Close
- Select Processing then Regression
 - o Select: Curve to be processed: v(t) ▼


Processing interval	Start: 0
	End: 2
Equation	V(t)
 - o Click on Plot

We read for example:

$23.9E-3 t + 182E-3$


Correlation coefficient:

0.2%

Line voltages

2.6 Set-up


2.7 Experiment

- Perform the assembly.
- Select the Test and calibrate icon
 - o Select: Channel: **Voltmeter /**
 Rated capacity: **⊙ 25 V**
⊙ Instantaneous

- o Adjust the zero
- o Click on OK
- o Select: Channel: Voltmeter//
Rated capacity: ☐ 25 V
☐ Instantaneous
- o Adjust the zero
- o Click on OK
- o Select: Channel: Voltmeter///
Rated capacity: ☐ 25 V
☐ Instantaneous
- o Adjust the zero
- o Click on OK
- Select the Parameters icon
 - o Validate: Acquisition mode, f(t) ▼
Acquisition time: 2 s ▼
Number of points: 501 ▼
☒ Synchronisation
 - o Select: Channels on the ordinate: Voltmeter/
Voltmeter//
Voltmeter///
 - o Click on Synchro...
 - o Select: Channel: Voltmeter/ ▼
Level: 0
Synchronisation ☒
☐ Increasing
 - o Click on OK
 - o Click on OK
- Select the Launch icon
 - o Select: Magnitude: u12
Names of the curves to be acquired: u12(t)
Magnitude: u23
Names of the curves to be acquired: u23(t)
Magnitude: u31
Names of the curves to be acquired: u31(t)

- o Click on Acquisition

2.8 Processing the curves


- Click on the Representation of curves icon
 - o Select: Scales: ☒ Automatic.
 - o Click on the Colours & Styles icon to modify the colour and the lines of the curves.

Represent $u_{12}(\pi)$, $u_{23}(\pi)$ and $u_{31}(\pi)$.

- Select Processing then Mapping
 - o Select: Curve to be mapped: u12(t) ▼

Abscissa	S.I. Units:	pi
	Expression:	2*t
Ordinate:	Name:	u12
	S.I. Units:	V
	Expression	u12
	Mapped curve:	u12(phi)
 - o Click on Plot then Save and then Close

Do the same with the other curves $u_{23}(t)$ and $u_{31}(t)$.


- Click on the Pointer icon
 - o Use the cross hair to measure:
- The period of each of the line voltages
- The offset from one line voltage to the other
- The maximum value of each of the line voltages.

Use Zoom for a more accurate reading.

Phase to ground and line voltages

2.9 Set-up


2.10 Experiment

- Perform the assembly.
- Start the computer; click on the Générís icon (the indicator lamp on the console comes on).
- Select the Test and calibrate icon
 - o Select: Channel: **Voltmeter /**
 Rated capacity: ☒ 25 V
☒ Instantaneous
 - o Adjust the zero
 - o Click on OK
 - o Select: Channel: **Voltmeter//**
 Rated capacity: ☒ 25 V
☒ Instantaneous
 - o Adjust the zero
 - o Click on OK
 - o Select: Channel: **Voltmeter///**
 Rated capacity: ☒ 25 V
☒ Instantaneous
 - o Adjust the zero
 - o Click on OK
- Select the Parameters icon
 - o Validate Acquisition mode: f(t) ▼
 Acquisition time: 2 ms ▼
 Number of points: 501 ▼
☒ Synchronisation
 - o Select: Channels on the ordinate: **Voltmeter/**
Voltmeter//
Voltmeter///


- o Click on Synchro...
 - o Select: Channel: **Voltmeter/** ▼
Level: 0
Synchronisation ☒ ☒ Increasing
 - o Click on OK
 - o Click on OK
- Select the Launch icon
 - o Select: Magnitude: v1
Names of the curves to be acquired: v1(t)
Magnitude: v2
Names of the curves to be acquired: v2(t)
Magnitude: u12
Names of the curves to be acquired: u12(t)
 - o Click on Acquisition

2.11 Processing the curves

- Click on the Representation of curves icon
 - o Select: Scales: ☒ Automatic.
- Click on the Colours & Styles icon to modify the colour and the lines of the curves.
- Click on the Pointer icon
 - o Use the cross hair to measure:
 - The period for each of the phase to ground and line voltages
 - The phase shift in line voltages u_{12} in relation to phase to ground voltage v_1
 - The maximum value of each of the voltages.

Use Zoom for a more accurate reading.

- Select Processing then Formula
 - o Select: Curve to be processed: v1(t) ▼
 - Ordinate: Name: U12
 - S.I. Units: V
 - Formula U = : v1-v2
 - Resulting curve: U12(t)

- o Click on Plot then on Save

Compare $u_{12}(t)$ and $U_{12}(t)$.

3 Viewing the voltages on an oscilloscope

3.1 Objective (see 2)


3.2 Equipment

- An extra-low voltage three-phase generator
- A dual-trace oscilloscope
- Four connecting wires.

3.3 Experiment


Work with a 150 Hz frequency to obtain a stable signal.

- ❖ View v_1 and v_2

- Press the INVERT Y2 button to view $-v_2$
- Press the ADD button to view $v_1 - v_2$
- Plot the oscillogram

- ❖ View u_{12}

Compare the two signals $v_1 - v_2$ and u_{12} .

To view the three voltages, use the two-way switch for phases 2 and 3.

4 Fresnel representation of voltages

We associate to every sinusoidal signal $u(t)$ a Fresnel vector $\vec{U}$ that rotates at an angular velocity ω .


Represent the Fresnel vectors $\vec{V}_1$, $\vec{V}_2$ and $\vec{V}_3$ associated to the three phase to ground voltages v_{1N} , v_{2N} and v_{3N} .

Represent the Fresnel vector $\vec{U}_{12}$ associated to the line voltage u_{12} knowing that: $u_{12} = v_{1N} + v_{N2}$

That is $\vec{U}_{12} = \vec{V}_1 - \vec{V}_2$

Graphically check the relation between the phase to ground and line voltage.

Such a system is called balanced.

5 After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your
disposal from Monday
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries