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***Material properties***

**Ref :  
242 011**


**Français – p 1**

**English – p 4**

**Version : 6010**

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***Elastic membrane manoscope***

## 1 Principe - description

### 1.1 Principe

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
## 1 Principal – description

### 1.1 Principal

The device enables the study of pressure distribution that exists in an liquid.

### 1.2 Description

The apparatus is made up of :

- a cylindrical capsule  $\varnothing$  50 mm for putting into the liquid. The pressure is exerted on a very supple rubber membrane covering the capsule and maintained by a moveable supple ring (11).
- a replacement membrane.
- a system of two co-axial tubes made from plastic :
  - the exterior tube in the shape of a J , linked to a panel supporting the manometer by a metal rod (1) with an open plastic ring (9), can be mounted or lowered into the liquid. It has silk screen printed markings every 50 mm (10).
  - the interior tube linking the capsule and the tube (7) open on the manometer can be turned using a knurled knob (8) : the movement is transmitted to the capsule by a waterproof supple connection (6).
- a manometer made up of a transparent plastic tube (4) mounted on a plate (2) with a scale in millimetres for using during experiments and a scale in centimetres which is particularly visible during classes (3).

- a 30 mL flask of coloured solution.

## 2 Set up

### 2.1 Fixing

Mount the apparatus vertically by blocking the metallic metal rod using a tightening nut on a base with rod (not supplied).

### 2.2 Refilling the manometer

1) Isolate the manometer by dividing the two supply plastic tubes. Pour the liquid into the U tube with a fine pipette, syringe, mercury funnel, etc...until the two levels reach more or less the zero mark (excess liquid can be ejected by softly blowing on one of the exits).

It is useless refilling the manometer to the closest millimetre; the final adjustment by simply pulling the U tube upwards or downwards will have the effect of reducing or increasing the length of the loop.

2) Link up the rubber tube with the bent part of the manometer: this operation will make the levels of liquid uneven.

3) Lightly hitch up using your thumb the rubber ring found just above the knurled ring covering the J tube to uncover the small hole (7) found in the lower tube. The manometer immediately returns to zero. Leave the rubber to fall to block off the hole once again.

## 3 Use

### 3.1 Necessary material

<i>Designation</i>	<i>Part no.</i>
Support with stem	701293
Nut of tightening	703452
Beaker F.B 3 L	713262

To completely take advantage to the device, it is recommended to use a recipient with a depth of at least 200 mm, so as to use the complete height of the J tube.

### 3.2 Experimentation

Sliding easily in the taut ring, the tube can be mounted or lowered, pivoting on itself in order to explore the complete water surface in the capsule.

The J tube can be stopped at a level, by turning the upper knurled ring, making the capsule pivot on a diametrical axis which shows that in these conditions, the mean pressure that it supports is independent to its position.

## 4 Upkeep

The device has practically no maintenance; none of the parts are sensitive to water. It is however advised to wash it after each use.

When not using the device, put talc on the rubber membrane and close the open manometer hole to stop water from evaporating.

The membrane can be changed after a certain amount of time (cracking, getting older). As this is foreseen, a replacement is supplied. To put it in place:

- Remove the plastic ring from around the capsule,
- Put back on the ring, the ring will be automatically taut again.

## 5 Replacement parts

<i>Designation</i>	<i>Part no.</i>
Membrane of replacement	242012
Fixing of the membrane	2 x 8574101

## 6 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

**JEULIN - TECHNICAL SUPPORT**  
**Rue Jacques Monod**  
**BP 1900**  
**27 019 EVREUX CEDEX FRANCE**  
**+33 (0)2 32 29 40 50**

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*Rubrique FAQ*

## Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

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