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Material properties

**Ref :
253 080**

Français – p 1


English – p 5

Version : 6010

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Thermal energy transfer kit

1 Description

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
1 Description

Each block has a hole into which the temperature sensor is inserted (not supplied).

This device is designed to study the transfer of thermal energy by conduction and radiation. It consists of a plastic kit that is used to store the equipment and for the experiment.

The kit contains the following elements placed in a preformed insulating foam:

- 2 parallelepipedal blocks in stainless steel of the same size
- 2 parallelepipedal blocks in aluminium also of the same size.
- 2 threaded shafts in plastic used for moving the blocks (especially when they are hot).

The kit may be completed by two Initio[®] thermometers, **part no. 251 040**,
a stopwatch **part no. 351 037**
and a water bath 20-80°C, **part no. 591 004**.

2 Technical characteristics

- Dimensions of the kit: 170 x 120 x 50 mm
- Dimensions of the metal blocks: 30 x 55 x 10 mm
- Ports for the temperature sensors: $\varnothing = 3$ or 6 mm

3 Objectives

Two types of heat transfers can be carried out:

- Transfer of thermal energy by conduction when the two blocks are in contact with one another
- Transfer of thermal energy by radiation (and a little by conduction due to the presence of air) when the two blocks are separated by a thin layer of air (about 2mm).

The blocks are placed in an adiabatic chamber (due to the powerful insulating properties of foam), their initial temperatures being different.

There is a transfer of thermal energy between them: the one with a higher temperature gives up part of its energy to the other. The temperatures of the two bodies tend towards a common limit, which is their equilibrium temperature.

4 Experiments

The experiment can be performed in two ways: manually or computerized.

4.1 Manual experiment

4.1.1 Implementation

One of the two metal blocks is heated by immersing it in very hot water at a temperature of 80°C for a few minutes so that the temperature of the block is uniform.

The thermometer probe must also be at 80°C to carry out correct initial measurements.

The other block is already placed in the slot provided in the foam of the kit at ambient temperature.

Transfer by conduction:

Place the two stainless steel blocks side by side in the left slot of the case. As they must be in proper contact with each other, there is some friction amongst them when they are inserted.

Transfer by radiation:

Place the two aluminium blocks in the right slot of the case. They are separated by 2 mm of air. Two matchsticks can be used to lock them in order to keep their separation constant.

In both cases, the blocks are covered with insulating covers pierced for the insertion of thermometers.

After a few seconds (to let the thermometers attain the temperature of the blocks), we note down the temperatures at constant intervals, for example 10 s for transfer by conduction and 1 min for transfer by radiation, which is slower.

Measurements of the two temperatures are taken simultaneously “on the fly” and then are entered in a table. The thermal equilibrium is attained after about 3 min for transfer by conduction and after 15 to 20 min for transfer by radiation.

4.1.2 Data and measurements


In partners, each student is responsible for noting down one of the temperatures; one of them reads the stopwatch and gives a signal every 10 s for the other to make a reading; the other is responsible for copying the values read in a table.

Here is an example of temperatures noted down manually in the case of transfer by conduction:

t (s)	0	10	20	30	40	50	60	70	80	90
θ_1 (°C)	79.0	74.5	69.6	66.1	63.6	61.1	59.2	57.6	56.3	55.2
θ_2 (°C)	26.0	33.0	36.5	39.1	41.2	43.0	44.4	45.5	46.4	47.1
t (s)	100	110	120	130	140	150	160	170	180	190
θ_1 (°C)	54.3	53.5	52.8	52.3	51.8	51.2	50.9	50.6	50.2	50
θ_2 (°C)	47.7	48.1	48.4	48.7	48.9	48.9	49.0	49.1	49.2	49.2

4.1.3 Processing of measurements

The transfer of experimental points and the plotting of temperature variation curves $\theta_1 = f(t)$ and $\theta_2 = g(t)$ can be represented either on a graph paper or using the Generis plotter as shown below.

The students can observe the existence of a common stagnation temperature, equilibrium temperature and note that this equilibrium is not established instantly. It can be estimated that around three minutes are required to attain it.


The shape of the curves in solid lines (result of a modelling) shows that the temperature variations are exponential. The lack of symmetry of the two curves in relation to the equilibrium temperature suggests that there is loss of heat (especially for the hotter solid).

The theoretical equilibrium temperature should have been the average of the initial temperatures (the two blocks are identical) i.e. 52.5 °C, and yet it is only at 49.5°C which is a relative deviation of about 6%.

4.2 Computerized experiment

4.2.1 Additional equipment required:

- 2 Initio[®] thermometers with analog outputs **part no. 251 045**
or 2 K type thermocouple probes **part no. 253 063**
+ 2 K type thermocouple adaptors **part no. 452 106**
- An ESAO[®] 4 (or 3) console **part no. 451 101**
+ an ESAO[®] 4 (or 3) card **part no. 451 300**
or an ESAO[®] 4+ console **part no. 451 401**


4.2.2 Implementation

The experimental conditions are identical to those of the manual experiment, only the entry of measurements is automatic.

It is started by the operator and programmed for example in the parameters of Generis for duration of 15 min, for example, for transfer by radiation.

The results obtained appear as in the reproduction of the screen below.

The red lines are the data points entered and the blue ones represent the exponential models. The equilibrium temperature is a line with 44.5 °C as ordinate.

Here too there is loss of heat: The two curves representing the variations in temperature are not symmetrical in relation to the θ_{equi} line and the equilibrium temperature must be the average of the initial temperature, i.e. 50°C because the blocks are identical. This time the relative deviation is 11%.

5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service * :
+ 33 (0)2 32 29 40 50

** Call cost not included*